

en

Folding Wheelchair

Model 9.050, Budget

Model 3.940, Format

Operating manual



ORTOPEDIA
A BRAND FROM MEYRA

Contents


Meaning of the applied markers	6
Introduction	6
List of models	6
Indications / contraindications	6
Acceptance	7
Intended purpose	7
Use	7
Wheelchair folding/unfolding/carrying	8
Adjustment	8
Combination with manufacturer foreign products	9
Reinstallment	9
Life span	9
Overview	10
Brake	11
Pressure brake - user	11
Locking the brakes	11
Releasing the brakes	11
Service brake	11
Drum brake for accompanying persons	12
Function as operating brakes	12
Leg supports	13
Calf belt	13
Removing the calf belt	13
Attaching the calf belt	13
Length adjustment of the calf belt	13
Lower leg support	14
Footplates	14
Leg support upper part	15
Turning the leg supports to the side	15
Swivelling in the leg supports	16
Removing the leg supports	17
Attaching the leg supports	17

Arm supports	18
Swivelling up the arm support	19
Swivelling the arm support inward	19
Swivelling down the arm support	19
Removing the arm support	20
Inserting the arm support	21
Back support	22
Adjustable back	22
Fitting the back belt	22
Wheels	23
Drive wheels	23
Support castors	24
Stick-in support castors	24
Removing/inserting the support castors	24
Support castor length	24
Retaining strap	25
Basic safety information	26
Accompanying person	26
Transfer out of the wheelchair	26
Reaching for objects	27
Driving on falling, rising or transverse gradients	27
Crossing obstacles	28
Loading the wheelchair	28
Transport in vehicles	29
Transport safety of the empty wheelchair	29
Transport of people inside a motor vehicle	29
Transport in public methods of transportation	29
Driving on public highways	30
Cleaning	30
Finish	31
Disinfection	31

Repairs	31
Repairs	31
Customer Service	31
Spare parts	31
Information for extended pauses of use	32
Disposal	32
Maintenance	32
Maintenance	32
Maintenance schedule	33
Technical data	35
Tyre pressure of pneumatic tyres	35
Applied norms	35
Data according to ISO for model 9.050 Budget	36
Further technical data for model 9.050 Budget	37
Data according to ISO for model 3.940 Format	38
Further technical data for model 3.940 Format	39
Meaning of the symbols on the washing instruction	40
Meaning of the labels on the wheelchair	41
Meaning of the symbols on the type plate	42
Inspection certificate	43
Warranty / Guarantee	44
Warranty / Guarantee section	45
Inspection certificate for transfer	45
Notes	46

MEANING OF THE APPLIED MARKERS

Safety instructions with a coloured background are mandatory and need to be observed under any circumstance!

-  This symbol indicates tips and recommendations.
- [] Reference to a picture number.
- () Reference to a function element within a picture.


INTRODUCTION

Read and observe this manual before first operation.

Children and juveniles should read this documentation together with their parents respectively a supervisor or accompanying person before first use.

Your specialist dealer should also introduce you into the use of your wheelchair.


This operating manual is to help you get accustomed to the handling of the wheelchair as well as to prevent accidents.

-  Please note that the illustrated equipment variants can deviate from your model.

We have therefore also listed chapters with options that might not be applicable for your vehicle. A list of the available options and accessories can be viewed in the order form of your wheelchair.

Users with visual impairments can find the PDF-files together with further information on our website:

< www.meyra.com >.

-  Contact your specialist dealer when required.

Alternatively users with visual impairments can have the documentation read out by a helper.

Inform yourself regularly about product safety and possible recalls of our products in the < *Infozentrum* > on our website: < www.meyra.com >.

We have developed a product that fulfils the technical and legal requirements of medical devices. For information about a severe accident that can still not be ruled out completely, please use our E-mail address < info@meyra.de > and inform the responsible governmental agency of your country.

LIST OF MODELS

This operating manual applies to the following models:

Model 9.050
Model 3.940

INDICATIONS / CONTRAINDICATIONS

In case of allergic reactions, skin rashes and/or pressure sores during the use of the wheelchair sores contact a doctor immediately.

In order to prevent contact allergies, we recommend to use the wheelchair only when wearing clothes.

The wheelchair serves as support in case of extremely pronounced walking disability resp. inability to walk, e.g. caused by:

- paralysis,
- defects or deformations of limbs,
- loss of limbs,
- damages or contractures in joints,
- other diseases (severe to fully pronounced impairments in mobility/

walking in cases of structural an/or functional damages of the lower extremities, u. o. amputations, results of injuries, musculoskeletal/neuromusculoskeletal caused movement disorders).

The wheelchair may not be used with a disability to sit.

The wheelchair may not be used without an accompanying person, in cases of:

- uncontrolled balance disorders,
- cognition disorders,
- loss of extremities on both arms,
- contractures / damages of the joints on both arms,
- insufficient or no eyesight.

🗨️ Please ask your doctor, therapist or specialist dealer to these and possible other risks in combination with this product.

ACCEPTANCE

All products are checked for faults in the factory and packed in special boxes.

- 🗨️ However, we request that you check the vehicle for possible transport damage immediately on receipt – preferably in the presence of the carrier.
- 🗨️ The packaging of the wheelchair should be stored for a further transport that might become necessary.

INTENDED PURPOSE

The wheelchair family serves for independent movement of one person or transport through an accompanying person.

USE

Do not reach into the spokes or grip the tyres of the rotating wheels. – Danger of injury!

Do not reach into the space between the pressure brake and the tyre when propelling the wheelchair. – Danger of injury!

Avoid a jerky propulsion of your wheelchair. – Danger of tipping over or tilting!

Do not use the wheelchair without attached leg supports and arm support units!

Your wheelchair with a seat width of more than 500 mm is not suited for transport in a train.

The wheelchair is applicable on level, firm surfaces and can be used as follows:

- for indoors (e.g. apartment, day care),
- outdoors (e.g. paved paths in parks),
- as a companion on tours (e.g. in a bus or train),
- 🗨️ In specific dimensions (overall height, overall width, overall length, turning radius) your wheelchair is not suited for transport within a train.
- 🗨️ Get information on possible limitations before the transport.
- transport within an airplane should only be done with PU-tyres.
- 🗨️ Ask the airline for possible limitations before transport in an airplane.
- Never expose the wheelchair to extreme temperatures and damaging environmental conditions, such as sunlight, extreme cold or salty water.
- Sand and other dirt particles can seize on moving parts and render them without function.

- Use only the handrims on the drive wheels to propel your wheelchair.

You must not let yourself be carried in your wheelchair by having the wheelchair lifted from the floor. Parts that are not securely fixed, e.g. side elements or leg supports, can come away and thus cause an accident.

Only apply the wheelchair within the scope of the specifications and limitation described in chapter *Technical data* on page 35.

WHEELCHAIR FOLDING/ UNFOLDING/CARRYING

Your wheelchair can be folded without the aid of tools.

1. Remove the seat cushion, if applicable.
2. Lift up the footplates or remove the leg supports.
3. Pull the seat belt upward in the centre of the front and back.
4. If necessary remove the connection rod of the double handrim drive or the push bar.

Your wheelchair can be carried without difficulty when folded.

For this, push one lower arm from the front under the upward folded seat belt. Place other hand underneath the rear seat fold for support.

Lift the wheelchair to horizontal position.

Tip the wheelchair slightly to one side in order to make the unfolding easier. At the side resting on the floor, press the seat tube down to the end stop.

- ☞ You may need to press down both seat tubes with the hands in order to do this.
- ☞ Ensure that the seat tube sits securely on the support surface.

ADJUSTMENT

Always have adaptation, adjustment or repair work carried out by a specialist dealer.

The wheelchair offers manifold adjustment possibilities to individual vital statistics. Before first use an adaptation of the wheelchair and a practical instruction in the functionalities of your wheelchair should be carried out by your specialist dealer. The adaptation will take into account the driving experience, the physical limits of the user and the main place of use of the wheelchair. Before first use, check the functionality of your wheelchair.

- ☞ We recommend a regular control of the wheelchair adjustment in order to ensure a long-term optimal provision even with changing illness/handicap patterns of the user. Especially for juveniles an adjustment every 6 months is recommendable.

COMBINATION WITH MANUFACTURER FOREIGN PRODUCTS

Any combination of your wheelchair with components not supplied by us generally results in an amendment to your wheelchair. Please inquire with us if there is a valid combination clearance/release.

Auxiliary drives

Mounting of additional- and front attachment drives constitute an extreme amendment of the wheelchair. The amendment will result in a forfeiture of the declaration of conformity for the wheelchair.

- ☞ The verification of conformity in combination with the valid medical device directive / medical device ordinance lies in the responsibility of the supplier of the additional or header drive. He is also responsible for the product included documentation of the medical product combination.

REINSTALLMENT

The wheelchair is suited for reinstallation. With the building block system the wheelchair can be fit to accommodate different handicaps body sizes. Before reimplementation the wheelchair is to undergo a complete inspection.

- ☞ Hygienical measures required for reinstallation are to be carried out according to a validated hygienic plan and must include disinfection.

The service manual, intended for the specialist dealer provides information to the reinstallation and reinstallation frequency of your wheelchair.


LIFE SPAN

We expect an average life span of about 4 years for this wheelchair, as far as the wheelchair is applied for its designated purpose and all maintenance and service guidelines. The life span of your wheelchair depends upon the frequency of use, the application environment and care. The implementation of spare parts can prolong the life span of the wheelchair. As a rule spare parts are available up to 5 years after production is discontinued.

- ☞ The indicated lifespan does not constitute additional guarantee.

OVERVIEW

The overview shows, representative for all models, the most important components of the wheelchair.



Pos. Description

- | | |
|---------------------------------|--|
| (1) Push handle | (8) Steering wheel |
| (2) Back support | (9) Locking lever- Arm support |
| (3) Arm support | (10) Locking knob – Quick release axle |
| (4) Seat belt/seat cushion | (11) Brake lever – pressure brake |
| (5) Locking lever – Leg support | (12) Handrim |
| (6) Calf belt | (13) Driving wheel |
| (7) Footplate | |

BRAKE

Arrange an immediate repair of the brakes by your specialist workshop if the braking performance reduces.

Wheelchairs with PU-tyres when stored for longer periods of time, should be put away safely without activated pressure brakes, since otherwise permanent deformations to the tyre surface could occur.

By locking the brakes with the brake lever (1), the wheelchair is secured against rolling away unintentionally (parking brake).

Depending on the version, the wheelchair can be equipped with pressure brakes (2) or with drum brakes (3).

Pressure brake - user

Locking the brakes

To secure the wheelchair against any unintentional rolling, press both brake levers forward all the way [4].

- It should not be possible to push the wheelchair forward when both brakes are locked.

Releasing the brakes

Pull both brake levers back all the way [2].

Service brake

The wheelchair is braked down with help of the handrims.

- If needed use suitable gloves in order to brake down the wheelchair.




Drum brake for accompanying persons

The drum brake is activated by the accompanying person through the brake levers (1).

Function as operating brakes

Use both brake levers evenly and only lightly in order to achieve a controlled deceleration of the wheelchair.




Locking the drum brakes

Pull both brake levers (1) [2] and press the locking latch (3) forward.

Release the brake lever. – Depending on the adjustment the latch engages in the first (A) or second (B) locking position and activates the brake.

⚠ If the third locking position (C) is required to activate the brake, it must be readjusted by a specialist workshop.



Loosen the drum brakes

Pull both brake levers (1) [2] until the latches (3) automatically jump out of the locking points.

Let go of both brake levers. – The parking brakes are released and the wheelchair ready for use.

LEG SUPPORTS

Before any actions on the leg supports the wheelchair is to be secured against unintentional rolling motions.

- Therefore observe chapter *Brake* on page 11.

Calf belt

Do not drive without the calf belt. – Danger of accidents.

The removable calf belt (1) prevents the feet from sliding off the back of the footboard.

- The calf belt must be removed in order to swivel away the leg supports.
- The calf belt is omitted for height adjustable leg supports and is replaced by a calf pad.

Removing the calf belt

For removal the calf belt is to be pulled from the attachment pins [2].

Attaching the calf belt

For attachment both loops of the calf belt are slid over the attachment pins [3].

Length adjustment of the calf belt

For length adjustment, the calf belt is guided around the special attachment pins [3] and adjusted in length with a velcro fastener.



Lower leg support

The footplates are to be folded up for entry into, exiting the wheelchair or scuttling (moving the wheelchair with help of the feet) [1].

- ☞ Check the locking points!
 - Remove both feet from the footplates.
 - Remove lower calf belt, if present.
 - ☞ Therefore observe chapter *Calf belt* on page 13.
- ☞ Before starting to drive the footplates are to be folded down again [2] and the calf belt attached.

Footplates

The footplates can be folded outward and up [2] resp. inward and down [3].



Leg support upper part

The upper leg support with an inserted lower leg support is termed leg support.

Turning the leg supports to the side

Leg supports turned to the side are released automatically and can easily come off. Note this when handling (e.g. transport).

For easy transfer out of/into the wheelchair as well as driving closer to a closet, bed or bathtub the leg supports can be swivelled away toward the in-/outside [1] and [2].

- ☞ Remove the calf belt, if applicable, before swivelling away the leg supports.
 - ☞ Therefore observe chapter *Calf belt* on page 13.

Fold up the footplates resp. footboard in order to swivel away the leg supports.

- ☞ Therefore observe chapter *Lower leg support* on page 14.

Afterwards pull or press the respective locking lever (3) backward and swivel the corresponding leg support outward [1]+[2].



Swivelling in the leg supports

For inward swivelling, let the leg supports swivel forward until the lock audibly engages [1].

- After audibly swivelling the leg supports inward check the respective locking device.
- Afterwards observe chapter *Lower leg support* on page 14.



Removing the leg supports

With removed leg supports there is an increased danger of tilting backwards.

For easy transfer into and out of the wheelchair as well as a reduced wheelchair length (important for transport) the leg supports can be removed [1].

- Remove the calf belt before swivelling away the leg supports.
 - Therefore observe chapter *Calf belt* on page 13.

For removal first swivel the leg support sideways [2] and then remove them toward the top [1].

- Therefore observe chapter *Turning the leg supports to the side* on page 15.
- Watch for possible danger of jamming!

Attaching the leg supports

For inserting press the leg supports, swivelled to the side, parallel to the front frame tube and lower it into place. – In doing so the holding pin must slide into the frame tube.

- After attachment swivel the leg supports inward [3].
 - Therefore observe chapter *Swivelling in the leg supports* on page 16.



ARM SUPPORTS

Do not use the arm supports to lift or carry the wheelchair.

Do not drive without the arm supports.

No not grab between the frame and arm support. – Danger of squashing!

When the wheelchair is being pushed by an attendant the user is to place his hands onto the arm cushions or in his lap and not at the sides between body and arm support. – Danger of squashing the fingers!

During lateral transfer without arm support, there is a danger of unintentionally releasing the brake and soiling the clothes on the drive wheel.

The arm supports [1] are removable [2] and simultaneously serve as padded arm support, clothes guard and wind protector.



Swivelling up the arm support

For transfer out of/into the wheelchair the arm support can be swivelled toward the back [1] as well as turned behind the back support [2].

Fold the front locking lever up in order to swivel up the arm support (3).

Afterwards swivel the arm support toward the back and up [1].



Swivelling the arm support inward

For inward swivelling, first turn the arm support upward [1] and then turn it toward the rear behind the back support [2].


- ☞ View chapter *Swivelling up the arm support* on page 19.
- ☞ On smaller back belt heights the arm support cannot be turned behind the back support.
- ☞ On the version with drum brake, the arm support cannot be folded behind the back support.



Swivelling down the arm support

In order to swivel down the arm support, carry out the sequence of chapters *Swivelling the arm support inward* on page 19 and *Swivelling up the arm support* on page 19 in reverse order.

- ☞ Afterwards check the locking device of the arm support and secure it!



Removing the arm support

Version 1:

To remove the arm support [1] fold the front and rear locking levers up (2).


Then remove the arm support in an upwards direction [1]

Version 2:

To remove the arm support [1] swivel the front locking lever up (2), then swivel the arm support up and inward.

☞ View chapter *Swivelling the arm support inward* on page 19.

Afterwards lift the arm support upward and off [1].






Inserting the arm support

First insert the arm support from the top into the rear bracket [1].

Then swivel the arm support down and press it into the front bracket [2].

- ☞ In doing so the rear tube of the arm support must lie in the guide groove on the back (3).

For safety fold the front and rear locking lever down (4)+(5).



BACK SUPPORT

Adjustable back

The tensioning of the back support [1] is adjustable.

Fitting the back belt

The adjustable back support strap has influence on the stability of the torso.


The pressure of the back must be spread evenly throughout the back cover.

The overlapping of the loop and hook strap has to be at least 10 cm.

To adapt the adjustable back strap, first pull the flap of the back support off and fold it forward [2].

Then open the Velcro fastener of the belt (3) that you wish to adjust and close it again after adjustment.

Afterwards fold the flap of the back support back down and fasten it with the velcro fastener [4].



WHEELS

Drive wheels

Do not use the wheelchair without mounted drive wheels.

After inserting the drive wheel the locking button (2) must stick a couple of millimetres out of the wheel nut.

The drive wheels are supported on a quick release axle [1] and can be removed resp. attached without tools.

First press the locking button (2) of the quick release axle in the center of the hub.

Afterwards remove or attach the drive wheel.

- ☞ The air pressure value for the tyres of the wheelchair can be read in the *Technical data* on page 35 or details on both sides of the tyre cover.
- ☞ If the drive wheels has too much side-ward lag or the quick release axle does not engage, contact your specialist dealer immediately for repair.
- ☞ No person may be seated in the wheelchair during assembly or removal. The wheelchair should stand on a level and firm surface. Before starting the disassembly work, support the frame to prevent the wheelchair from tipping over and secure it to prevent an unwanted movement or tipping over.




SUPPORT CASTORS

Support castors do not provide sufficient protection against tipping over in certain situations.

Definitely refrain from the following:

- Leaning the upper body far back.
- Jerky acceleration, especially when driving uphill.

In order to raise the tilting stability never drive without the support castor swivelled to the back.



A bent tube with two small wheels serves to increase tilting stability [1].

Stick-in support castors

The insertable support castors can be inserted from the back into the lower frame tubes [1].

Removing/inserting the support castors

Depress the spring button (3) in order to remove/insert a support castor.

Support castor length

In order to adjust the length of the support castors, press down the spring button (2) and telescopically pull or push the support castor to the desired length.

RETAINING STRAP

Make sure that no objects are trapped between belt and the body!

The retaining strap is not part of the retaining system for the wheelchair and/or the driver during transport in motor vehicles.

The retaining strap [1] is screwed from the back onto the respective back support tube. The retaining strap serves to stabilise the sitting position and prevents falling forward out of the wheelchair.

To fasten the retaining strap, pull both ends forward and audibly let the buckle click into place. To open the retaining strap, press down the release button and pull the two ends of the strap apart.

- ✎ The retaining straps can be adjusted in length and should not be pulled too tight.



BASIC SAFETY INFORMATION

This safety information is an extract of the *Safety and general handling instructions*, that can be found on our website: < www.meyra.com >.

Do not insert fingers into open frame tubes (for example after removing the arm supports, leg supports or support castors). – Danger of injury!

Scuttling (rolling motion of the wheelchair with help of the feet) is only permitted with creeping rolling speed on horizontal, level surfaces indoors.

When scuttling watch for the possibility of jamming the feet.

In doing so the sitting position may not be slid onto the front third of the seat.

A stable sitting position is to be kept while using the wheelchair, even when not in motion and especially on hills and slopes. – Danger of accidents!

In a safe sitting position the back of the user lies directly on the back support upholstery and the hip of the user is at the back end of the seat.

Transit out of the wheelchair on hills/slopes may only be carried out in emergencies and with the aid of an accompanying person and/or helper! – Danger of accidents!

Adjust the seat inclination only when the wheelchair is standing on a horizontal, level surface. A danger of tipping over exists on gradients!

Increased danger of tipping over when using the angle adjustable back support.

Before starting to drive make sure that no negative seat inclination has been adjusted, respectively that a safe seating position is ensured.

You should not smoke while using the wheelchair.

Exposure to direct sunlight can cause seat covers/upholstery, arm support pads, leg supports and handles to heat up to over 41 °C. – Contact with exposed skin can result in injury! Prevent such heating by parking the wheelchair in a shaded area.

Special attachment points for fastening objects that are to be carried along are not provided.

Accompanying person

The accompanying person must be made aware of all possible danger situation before the start of his/her supportive involvement. The parts of your wheelchair that are held onto by the accompanying person are to be checked for tight fit.

Transfer out of the wheelchair

Drive with the wheelchair as closely as possible to the spot where you want to switch out of the wheelchair.

- ☞ Herefore additionally observe chapters *Brake* on page 11 and *Swivelling up the arm support* on page 19.
- ☞ We recommend to conduct the transfer from the wheelchair together with an aid.

Reaching for objects

Avoid an extreme forward or backward inclination of the upper body when picking up or placing heavy objects. – Danger of wheelchair tipping over, especially in the case of narrow seat widths and high seat heights (seat cushion)!

Driving on falling, rising or transverse gradients

For safety reasons, the maximum permitted gradient is limited because the tip-over stability and the braking and steering behaviour are impaired by a reduced floor/road.

☞ Observe chapter *Technical data* on page 35.

For your own safety, we recommend that you use the assistance of an accompanying person for driving slopes of more than 1% tilt.

☞ This low gradient value is valid for persons that are not able to effect even the smallest centre of gravity change by changing the position of their upper body.

The accompanying person should be aware that the control of the wheelchair requires more strength on a downward slope/ hills than on a level driving surface and that there is an increased danger of slipping.

Never lean towards the downhill direction when driving on rising, falling or transverse gradients.

Avoid jerky changes of the driving condition (especially with critically adjusted driving parameters as for example an active wheel position).

Never drive faster than walking speed.

There is an increased danger of tilting while driving backwards.

While driving in curves and when turning on inclinations and slopes there is a danger of tilting.

The braking force transferred to the driving surface is much less on a downward slope than on a level driving surface and is further reduced by poor road conditions (e.g. rain, snow, grit, dirt). A dangerous slipping of the wheels due to excessive braking and an associated unwanted course deviation must be avoided by way of a careful dosed braking.

At the end of the downward slope, take care that the foot plates do not make contact with the ground and endanger you through a sudden braking effect.

Transverse surfaces to the driving direction (e.g. transversely sloped pavements) effect a turning of your wheelchair in the downhill direction. You or an accompanying person must compensate for this drift by a counter-steering.

Crossing obstacles

After having overcome stairs, support castors that were previously removed need to be mounted again resp. support castors that were swivelled to the side, swivelled back into place.

Each crossing of obstacles involves a risk! Even flat obstacles (e.g. sills, carpet edges, small steps or similar) represent a danger of tipping.

The crossing of obstacles is a special danger situation in which a combination of the safety advice in the sections headed uphill driving, downhill driving and driving transverse to a slope must be observed in addition to other safety advice.

For your own safety all obstacles should only be driven over with an accompanying person. The pedal tube is to be used as support for each manual tilting of the wheelchair.

The unassisted, safe, crossing of small obstacles (< 4 cm) is only possible with a sufficiently safe control of the wheelchair.

Always approach small obstacles, e.g. steps or curbs, slowly and at a right angle (90°) until the steering wheels almost touch the obstacle. Briefly stop the wheelchair and then drive over the obstacle.

Always drive over an obstacle forwards (not backwards) and simultaneously with both front/rear wheels in order to prevent a sideways tipping of the wheelchair.

Always maintain a safety distance between the wheelchair and drops, stairs and similar obstacles sufficient for reaction, braking and turning.

If possible, let one or more helpers lift you out of the wheelchair and carry you to the destination point.

The crossing of stairs is only possible with the assistance of two helpers who are familiar with the danger situation and firmly guide the wheelchair.

☞ The helpers may only hold (without lifting) the wheelchair with the push handles and the side frame.

Existing support castors are to be removed, since they might pose a danger to the aid under certain circumstances.

Loading the wheelchair

Do not use the leg supports, arm supports or accessory parts in order to lift the wheelchair!

Before lifting, the wheelchair is to be secured against unintentional rolling motions!

☞ Therefore observe chapter *Turning the leg supports to the side* on page 15.

The parts detached for loading must be carefully stowed and carefully attached again before the next journey.

The weight of the wheelchair is reduced when you remove detachable components.

The wheelchair can be also loaded with the aid of ramps or lifting platforms.

☞ Therefore observe the respective user manual and the load capacity of the ramp resp. lifting platform.

Transport in vehicles

The parts detached for the transport must be carefully stowed and carefully attached again before the next journey.

No special carrying points are allocated for carrying detachable components.

The following procedures may be necessary due to lack of space for the transport in vehicles:

- Remove the leg supports.
- Remove the arm supports.
- Remove or fold over the back support.
- Detach the drive wheels.

Transport safety of the empty wheelchair

All regulations and directions of the respective transport company are to be observed. – Ask for these before the transport.

The wheelchair is to be stored and secured in such a fashion that it does not present any danger to you or others.

Transport of people inside a motor vehicle

Whether your wheelchair is approved as a seat in a car, can be determined in chapter *Meaning of the labels on the wheelchair* on page 41.

- ☞ Observe the guideline < *Safety with Meyra-wheelchairs, also during transport in motor vehicles* >! – This document and further information are available in the < *Information center* > on our website < www.meyra.com >.

Transport in public methods of transportation

Your wheelchair is not designed for user transport in public transportation vehicles. Limitations may occur. We recommend use of one of the firmly built in seats of the public vehicle.

Should it nevertheless become necessary to carry out the transport while sitting in the wheelchair, the following needs to be observed:

- Use the space designated by the public transportation services for parking.
- Observe the regulations of the transport company before parking the wheelchair.
- Park your wheelchair opposite to the driving direction in the reserved space.
- The wheelchair is to be placed so that the back support it will be supported by the border of the parking space.
- One side of the wheelchair must also lie against the border of the parking space, so that the wheelchair cannot slide away in case of an accident or sudden braking manoeuvre.
- Additionally activate the parking brakes.

Driving on public highways

Observe the valid regulations for public traffic of your country and if necessary ask your specialist dealer for required accessories.

Your wheelchair can optionally be fit with lighting equipment. The lighting equipment consists of:

- Reflectors on the back support
- Reflectors on the drive wheels.

With limited visibility and especially in the dark we recommend to mount active lighting equipment and to turn it on in order to see better and be seen.

- ☞ When participating in public traffic the user is responsible for the functional and operationally safe condition of the wheelchair.
- ☞ The valid traffic regulations must be observed and abided when participating in public traffic.
- ☞ Wear light-coloured and conspicuous clothing when driving in darkness.
- ☞ When driving in the dark avoid using the road or bicycle lanes.
- ☞ Observe that the lighting equipment is not covered by clothes or any other objects attached to the wheelchair.
- ☞ In case of physical handicaps, such as blindness the use of the wheelchair is only permitted with an accompanying person.

CLEANING

The plastic panelling is attacked through non-ionic tensides as well as solvents and especially alcohol.

The cushions and covers are normally fit with care instructions (instruction for care).

- ☞ For this observe chapter *Meaning of the symbols on the washing instruction* on page 40.

In all other cases the following information is true:

- ☞ Clean the upholstery with warm water and hand washing liquid.
- ☞ Remove spots with a sponge or a soft brush.
- ☞ Wash off persistent dirt with commercial fine detergent.
- ☞ Do not soak! Do not machine wash!

Follow-up with clean water and allow to dry.

The chassis, the lighting equipment and wheels can be cleaned damp with a mild detergent. Afterwards dry off well.

- ☞ Check the chassis for corrosion damages as well as other damages.
- ☞ Only clean the plastic parts with warm water and neutral detergent or soft soap.

Further information to cleaning can be found in the < *Infozentrum* > on our website: < www.meyra.com >.

Finish

The high quality finish ensures an optimum of protection against corrosion.

- ☞ Should the coating be damaged with scratches or similar, these areas can be touched up with our paint pen available at the specialist dealer.

Slight lubrication of moving parts will ensure for their long functioning.

Disinfection

If the product is used by more than one person (for example in a care centre), the use of a commercial disinfectant is mandatory.

- ☞ Before disinfection the upholstery and handles are to be cleaned.
- ☞ A spray- or wiping disinfection is permitted with tested and accredited disinfectants.

You can get information on tested and permitted disinfectants and procedures at your national facility for health protection.

- ☞ During the use of disinfectants it can happen that surfaces might be affected in such a fashion that the long term functionality of parts can be limited.
- ☞ In doing so the manufacturers instructions are to be observed.

REPAIRS

Repairs are generally to be carried out by a specialist dealer.

Repairs

Trustingly contact your specialist dealer for maintenance work. He has been introduced to the maintenance.

Customer Service

In case you have any questions or need help please contact your specialist dealer who can assume counselling, customer service and repairs.

Spare parts

Spare parts can only be ordered from specialist dealers. In case of repair work, only original spare parts are to be used!

- ☞ Spare parts from other manufacturers can cause malfunctions.

The spare parts list with the respective part numbers and drawings is available at the specialist dealer.

In order to ensure the correct delivery of a spare part, always quote the corresponding serial number (SN) of the wheelchair! You will find this on the type plate.

Whenever a wheelchair repair is carried out by the specialist dealer, the supplementary information, e.g. assembly/operating instructions must be attached to the operating manual for the wheelchair, the date of the repair must be recorded and stated when ordering spare parts.


This should prevent wrong order details on future spare parts orders.

Information for extended pauses of use

No special measures are required during extended pauses of use. Nevertheless storage temperatures are to be observed.

- ☞ For this observe chapter *Technical data* on page 35.

DISPOSAL



The disposal must comply with the respective national law.

Please enquire about local disposal arrangements at your municipal authority.

MAINTENANCE


An incorrect or neglected cleaning and maintenance of the wheelchair results in a limitation of the product liability.

Maintenance

The following maintenance Instruction gives you a guide for carrying out the maintenance work.

- ☞ The maintenance plan does not give information about the actual extent of work determined on the wheelchair.
- ☞ In all cases of inspection or maintenance there is a danger that you injure yourself or others unintentionally.
– Therefore all inspection and maintenance is to be carried out with the indicated care.

Maintenance schedule

WHEN	WHAT	REMARK
Before starting out	General Test for faultless operation.	Carry out test yourself or with a helper.
Before starting out	Test brakes for faultless operation. Activate brake lever to the limit.	Carry out test yourself or with a helper. The locked wheels should not be able to turn under operating conditions. If they can still turn, the brakes must be repaired by an authorised specialist workshop.
Before starting out	Check pressure brake for wear. Move brake lever to the side	Carry out tests yourself or have a helper do it. If you notice any increasing slackness on the brake lever take the wheelchair to your specialist workshop immediately for repairs. – Danger of accidents!
Before starting out (when applicable)	Check air pressure of the tyres. Tyre filling pressure:  View Tyre pressure of pneumatic tyres on page 35.	Carry out test yourself or with a helper. Use a tyre gauge.
Before starting out	Check tyre profile.	Carry out visual check yourself. If the tyre profile is worn down or if the tyre is damaged, consult a specialist workshop for repairs.

WHEN	WHAT	REMARK
Before starting out	Check the back tubes and frame tubes for damages.	Carry out the test yourself or by a helper. If deformations or cracks occur in the welding seams, contact a specialist workshop immediately for repairs. – Danger of accident!
Especially before driving in the dark (when applicable)	Check the lighting Check the light- and indicator signal equipment as well as reflectors for immaculate performance.	Carry out test yourself or with a helper.
Every 8 weeks (depending on frequency of use)	Lubricate the following components with a few drops of oil. – Moving parts of the locking mechanism. – Brake lever bearings.	Do it yourself or with the aid of a helper. Components must be free from used oil residues before lubrication. Please ensure that excess oil does not contaminate the environment (e.g. your clothing).
Every 8 weeks (depending on frequency of use)	Check all screw connections for secure fit	Carry out test yourself or with a helper.
Every 6 months (depending on frequency of use)	Check – Cleanness. – General condition.	See Care. See Repairs.
Manufacturer recommendation: Every 12 months (depending on frequency of use)	Safety inspection – Wheelchair.	To be carried out by the specialist dealer.

TECHNICAL DATA

All data given in the < *Technical data* > refers to the standard version.

The overall length depends on the position and size of the drive wheels.

If not noted otherwise the dimensions are determined with drive wheels of \varnothing 610 mm (24").

The widths were determined with a handrim distance of 15 mm.

Dimensional tolerance ± 15 mm, $\pm 2^\circ$.

Short form of wheelchair dimensions:

SH = Seat height

SW = Seat width

SD = Seat depth

BH = Back support height

Calculation of the max. user weight:

The maximum total load is calculated on the basis of the unloaded weight of the wheelchair and the maximum passenger weight.

Additional weight due to subsequent additions or luggage reduce the maximum permissible passenger weight.

Example:

A driver wishes to take luggage with a weight of 5 kg. Thus, the maximum user weight is reduced by 5 kg.

Tyre pressure of pneumatic tyres

Maximum tyre pressure is printed on the tyres on each side.

Full tyre pressure – steering wheel

Standard:

2.5 - 3.5 bar = 36 - 50 psi

Full tyre pressure – drive wheel

Standard:

3.0 - 4.0 bar = 44 - 58 psi

Ultra-light running tyres:

6 bar = 87 psi

High pressure:

8 bar = 116 psi

Applied norms

The wheelchair complies with the norm:

- EN 12183: 2014
 - ☞ The weight of the dummy was 130 kg.
- ISO 7176-8: 2014
- ISO 7176 -19: 2008
 - ☞ Assessment of the Crashtest, in which the wheelchair is attached to the retaining system of the vehicle, has been carried out according to the testing methods of annex D.

The applied parts and components we use are in compliance to EN 1021-2 for resistance against inflammation.

Data according to ISO for model 9.050 Budget

	min.	max.
Overall length with leg support	1020 mm	1080 mm
Overall width	560 mm	660 mm
Overall dimensions	– kg	149,6 kg
User weight (incl. additional load)	– kg	130 kg
Weight of the heaviest part	– kg	11 kg
Actual seat depth	– mm	430 mm
Actual seat width	380 mm	510 mm
Folding length with leg supports	1020 mm	1080 mm
Folding width	– mm	280 mm
Folding height	960 mm	960 mm
Seat surface height at front edge (without cushion)	470 mm	510 mm
Seat angle	0°	4°
Back support angle	–°	90°
Back support belt height	– mm	420 mm
Foot support to seat (lower shank length)	380 mm	500 mm
Static stability downhill	–°	6°
Static stability uphill	–°	6°
Static stability lateral	–°	6°
Angle leg support - seat surface	–°	111°
Arm support height from seat surface	– mm	230 mm
Back support to front edge of arm support	– mm	280 mm
Handrim diameter	– mm	535 mm
Axle horizontal position	-30 mm	30 mm
Obstacle height	30 mm	100 mm
Minimal turning radius (normative recommendation according to ISO 1000 mm)	1250 mm	– mm
Pivot area	1300 mm	– mm

Further technical data for model 9.050 Budget

	min.	max.
Length without leg support	770 mm	830 mm
Overall height	– mm	960 mm
Seat cushion thickness	30 mm	60 mm
Length without leg supports, drive wheels (Support castors are removed or swivelled underneath the seat)	– mm	750 mm
Additional load	– kg	10 kg
Empty weight	19,6 kg	– kg
Transport weight (without cushion, drive wheels)	9 kg	– kg
Weight seat cushion	0,9 kg	– kg
Weight drive wheels (1 piece)	– kg	2,5 kg
<u>Steering wheel</u>		
ø 180 x 40 mm (7")		PU (puncture safe)
<u>Driving wheel</u>		
ø 610 mm (24 x 1 3/8")		PU (puncture safe)
<u>Climatic data</u>		
Ambient temperature		-20 °C to +50 °C
Storage temperature		-30 °C to +60 °C

Data according to ISO for model 3.940 Format

	min.	max.
Overall length with leg support	1020 mm	1080 mm
Overall width	560 mm	660 mm
Overall dimensions	– kg	147.5 kg
User weight (incl. additional load)	– kg	130 kg
Weight of the heaviest part	– kg	11 kg
Actual seat depth	400 mm	430 mm
Actual seat width	380 mm	510 mm
Folding length	1020 mm	1080 mm
Folding width	280 mm	– mm
Folding height	960 mm	– mm
Seat surface height at front edge (without cushion)	430 mm	510 mm
Seat angle	0°	4°
Back support angle	–°	90°
Back support belt height	400 mm	440 mm
Foot support to seat (lower shank length)	380 mm	500 mm
Static stability downhill	–°	6°
Static stability uphill	–°	6°
Static stability lateral	–°	6°
Angle leg support - seat surface	–°	111°
Arm support height from seat surface	– mm	230 mm
Back support to front edge of arm support	– mm	280 mm
Handrim diameter	485 mm	535 mm
Axle horizontal position	-30 mm	30 mm
Obstacle height	0 mm	100 mm
Minimal turning radius (normative recommendation according to ISO 1000 mm)	1250 mm	– mm
Pivot area	1300 mm	– mm

Further technical data for model 3.940 Format

	min.	max.
Length without leg support	770 mm	830 mm
Overall height	– mm	960 mm
Seat cushion thickness	30 mm	60 mm
Length without leg supports, drive wheels (Support castors are removed or swivelled underneath the seat)	– mm	750 mm
Additional load	– kg	10 kg
Empty weight	17,5 kg	– kg
Transport weight: (without leg supports, arm supports, cushion, drive wheels)	9 kg	– kg
Weight seat cushion	0,9 kg	– kg
Weight leg supports (1 piece)	– kg	2,1 kg
Weight drive wheels (1 piece)	– kg	2,5 kg
<u>Steering wheel</u>		
ø 142 x 37 mm (5,5")		Soft solid rubber
ø 180 x 40 mm (7")		PU (puncture safe)
<u>Driving wheel</u>		
ø 559 mm (22 x 1 3/8")		PU (puncture safe)
ø 610 mm (24 x 1 3/8")		PU (puncture safe)
<u>Climatic data</u>		
Ambient temperature		-20 °C to +50 °C
Storage temperature		-30 °C to +60 °C

Meaning of the symbols on the washing instruction

(the symbols correspond to European standard)



Wash as delicates with the indicated maximum temperature in °C.



Wash as regular laundry with the indicated maximum temperature in °C.



Do not bleach.



Not suited for the dryer.




Do not iron.




Do not dry-clean.

Meaning of the labels on the wheelchair



Attention!

Read the operating manuals and other provided documentation.



Do not lift the wheelchair at the arm supports or leg supports. Removable parts are not suitable for carrying.

Achtung
Bremsen nachstellen.
Ident.-Nr. 8390658

Attention

Readjust the brakes.

Achtung
Erhöhte Kippgefahr auf Steigung / Gefälle
besonders in Verbindung mit kurzem Rad-
stand.
Ident.-Nr. 205674400

Attention

Increased danger of tilting when on inclinations / slopes, especially in combination with short wheel base.



The product is **not** approved as a seat within a motor vehicle.



Indication for danger of jamming. – Do not reach in here.

Meaning of the symbols on the type plate



Manufacturer



Order number



Serial number




Production date



Permitted user weight



max. permissible total weight



Permitted axle weights



Max. permissible rising gradient



Max. permissible falling gradient

max. ... km/h

Permitted maximum speed



The product is approved as a seat within a motor vehicle.



Max. permitted user weight if the product is approved as a seat within a motor vehicle.



The product is **not** approved as a seat within a motor vehicle.

INSPECTION CERTIFICATE

Vehicle data:

Model:

Delivery note no.:

Serial-no.(SN):

Recommended safety inspection 1st year (at least every 12 months)

Stamp of specialist dealer:

Signature: _____

Place, date: _____

Next safety inspection in 12 months

Date: _____

Recommended safety inspection 2nd year (at least every 12 months)

Stamp of specialist dealer:

Signature: _____

Place, date: _____

Next safety inspection in 12 months

Date: _____

Recommended safety inspection 3rd year (at least every 12 months)

Stamp of specialist dealer:

Signature: _____

Place, date: _____

Next safety inspection in 12 months

Date: _____

Recommended safety inspection 4th year (at least every 12 months)

Stamp of specialist dealer:

Signature: _____

Place, date: _____

Next safety inspection in 12 months

Date: _____

Recommended safety inspection 5th year (at least every 12 months)

Stamp of specialist dealer:

Signature: _____

Place, date: _____

Next safety inspection in 12 months

Date: _____

WARRANTY / GUARANTEE

Failure to observe the instructions in the operating manual, improperly carried out maintenance work and, especially, technical changes and additions (add-ons) carried out without our prior consent will lead to a general loss of guarantee and product liability.

We accept legal liability for this product within the scope of our general terms and conditions and warranty and in certain cases other verbal resp. agreed upon guarantees. For warranty and guarantee demands please contact your specialist dealer with following Warranty/Guarantee section and the there included information on model description, delivery note number with delivery date and serial number (SN).

The serial number (SN) can be read off of the type plate.

Precondition for the acceptance of liability in any case is the intended use of the product, the use of original spare parts by authorised dealers as well as maintenance and inspections in regular intervals.

Guaranty is not granted for surface damages, tyres of the wheels, damages due to loosened screws or nuts as well as worn out attachment holes due to frequent assembly work.

Furthermore, damage to the drive and electronics caused by improper cleaning using steam cleaning equipment or the deliberate or accidental flooding of the components are also excluded.

Interferences through radiation sources such as mobile phones with high transmission power, HiFi-equipment and other extreme interference radiators outside of norm specifications cannot be declared as warranty or guarantee claims.

This operating manual as a part of the product is to be handed out in case of a change of owner.

For evaluation of our products you can use our < *Infozentrum* > sector < *PMS* > on our website < www.meyra.com >.

We reserve the right to make technical improvements.



This product is conform with the EC-Directive 93/42 EEC and fulfills the requirements of the from May 25th 2021 valid regulation (EU) 2017/745 for medical devices.

Warranty / Guarantee section

Please fill out! Copy if necessary and send the copy to the specialist dealer.

Warranty / Guarantee

Model designation:

Delivery note no.:

SN (view type plate):

Date of delivery:

Stamp of the specialist dealer:

Inspection certificate for transfer

Vehicle data:

Serial-no.(SN):

Model:

Delivery note no.:

Stamp of specialist dealer:

Signature: _____

Place, date: _____

Next safety inspection in 12 months

Date: _____

NOTES

Your specialist dealer

MEYRA GmbH

Meyra-Ring 2



32689 Kalletal Kalldorf
GERMANY



Tel +49 5733 922 - 311

Fax +49 5733 922 - 9311



info@meyra.de

www.meyra.de
