

Invacare® | **Flo-tech Lite Visco**



Low profile visco elastic low memory foam cushion offering exceptional comfort

The *Invacare Flo-tech Lite Visco* cushion is high density and temperature sensitive to the body, allowing it to conform to the user's shape. This affords exceptional levels of comfort whilst reducing the occurrence of peak pressures. This low profile cushion format is ideal for both chair and wheelchair use, facilitating side transfers.



Yes, you can.®

Low memory foam designed to maximise body contact area

The temperature-sensitive foam conforms to the user's body, maximising the contact area, delivering exceptional levels of comfort and optimising pressure reduction.

Contoured seating surface withommel and side support

Offers increased lateral support and pelvic stability to enhance the user's posture.



Features and Options



Ingress Resistant Cover

Cover is two-way stretch, water-resistant and vapour-permeable. Cover features an anti-slip base with carry handle.



Available in 61 sizes

The *Invacare Flo-tech Lite Visco* is available in a range of 61 sizes, from 25 x 25 to 48 x 53 cm (250 x 250 to 480 x 530 mm). (See size chart for details).



Three foam densities available

Available in three densities to accommodate a range of user weights – Soft (up to 32 kg), Medium (up to 51 kg) and Firm (up to 70 kg).

Technical data



A range of 61 sizes



Height: 7.8 cm (78 mm)



Weight: 1.26 kg



Max. User Weight: 100 kg



Cleaning: 80°

Flo-tech Lite Visco

Water Resistant
Vapour Permeable

Manufactured to European Standards:
Interior Foam: BS EN 597-1, BS EN 597-2, BS 5852:
1982 Crib 5
Cover: BS 7177: 1989 Crib 5

INVACARE Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5HZ
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
www.invacare.co.uk
Sales Order E-mail:
ordersuk@invacare.com

INVACARE Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2011 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare reserves the right to alter product specifications without prior consultation.

Flo-tech Lite Visco - EU - 06/2011

Colours



Black Smooth



Black Wipedown



Black Towelling



Yes, you can.®