


# Clou

9.500

 **Allowed user weight**  
120 kg





Long-life LED lighting available as an option



Folded



30° back adjustment can be selected as option



Adjustable seat and back

## The small folding one

### EFFICIENCY

- Home and inpatient application
- Great life span through precision steel tube frame and high tensility of important components
- Optimal reinstallment through simple adaptation possibilities to the user
- Different seat widths in one frame
- Long term safety in spare part provision through a wide spread specialist dealer network
- Little need for repairs by implementing high value and durable parts

### THERAPEUTIC BENEFIT

- Ideal for use in confined spaces due to compact dimensions
- Driving electronics adjustable in length and height, enabling comfortable, fatigue-free driving in an ergonomic posture
- Interactive programmable driving characteristics ensure individual adaptation
- Easy-to-operate control buttons enable handling even with very little strength
- Close driving to objects and easy transfer through leg supports removable to the side

### FIELD OF USE

- For users that are no longer able to move around with a manual wheelchair, even indoors
- For extreme functional limitations of the lower extremities and little hand and arm strength
- For driving in confined spaces
- As transport aid for operation through an accompanying person

ORDER

COST ESTIMATE

CUSTOM DESIGN

Please send the order form by e-mail to:

[info@meyragroup.com](mailto:info@meyragroup.com)

Customer service: +49 571 93292 555

Billing address / Customer ref. no.:

Delivery address:

Tax-ID

Consignment / Comments:

Order quantity:

**MEYRA**



## Allowed user weight

120 kg



## Hints

Our Terms and Conditions as stated at [www.meyra.de/AGB](http://www.meyra.de/AGB) shall apply

You can find information on the measurement process at [www.meyra.de/messverfahren](http://www.meyra.de/messverfahren)

Valid from 01.01.2025 – 31.12.2025.

Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at [www.kisss-by-meyra.de](http://www.kisss-by-meyra.de)

You avoid queries and incorrect orders thanks to the buildability check

## USER DATA

First name

Name

Street

ZIP CODE

City

Country

Birthday

Diagnosis 1

Diagnosis 2

Tel.

Mobil

E-mail

Order number

Health insurance fund / cost unit

Specialist trade Contact person/phone.

Supply specialist MEYRA

Approval info from week

## ANATOMICAL DATA (IN MM)

1) Seat depth

2) Seat width

3) Chest width

4) Chest depth

5) Elbow height

6) Armpit height

7) Back height


8) Head height

9) Forearm length

10) LLL

Body height

Weight (kg)



Remarks

I accept the Terms and Conditions of Meyra Distribution GmbH (valid as of January 2023)

<https://www.meyra.com/gtc/>

Date

City

Signature



## Technical specifications

<sup>1)</sup> With standard legrest CODE 93/808<sup>2)</sup> under test conditions with drive battery 38 Ah (5h), 45 Ah (20h)

Seat widths in mm	380 - 500 / 430 - 550
Seat depth in mm	430 <sup>1)</sup>
Back support height in mm	400
Lower shank length in mm	360 - 500 <sup>1)</sup>
Arm support height in mm	190 - 300
Front seat height in mm	510
Length with leg supports in mm	1100 <sup>1)</sup>
Length without leg supports in mm	810
Width min./max. in mm	600 / 700
Width folded in mm from	380
Height without head support in mm	910
Transport length in mm	810
Transport height min. in mm	910
Front wheel size in mm	200 x 50 (8")
Rear wheel size in mm	320 x 60 (12.5")
Max. obstacle height in mm	60
Turning radius in mm	1000
Transport empty weight basic configuration min.	min. 73 kg
User weight	max. 120 kg
Permitted overall weight incl. options	210 kg
Max. additional loading	10 kg
Driving speed	6 km/h
Motor power	2 x 180 W
Range	30 km
Permitted inclination/slopes	12%
Form of steering	indirect
Drive	Rear

☐ Please select by ticking the relevant boxes☒ Special construction/individual solution o.r. - on request☐ Standard version: If no other selection is made, the standard version will be delivered.**00 Chassis/Motorisation**

CODE	SPEED	MOTOR POWER	RANGE
115	6 km/h	2 x 180 W	max. 30 km

**01 Seat width (SW)**

CODE	SW	DESCRIPTION
38	SW 380 - 500 mm	Seat width, continuously adjustable through the arm supports from 380 - 500 mm, (seat surface width 430 mm)
43	SW 430 - 550 mm	Seat width, continuously adjustable through the arm supports from 430 - 550 mm, (seat surface width 430 mm)


**02 Seat system**

CODE	DESCRIPTION
738	Standard seat belt without cushion

**Seat system options**

CODE	DESCRIPTION
560	Seat cushion 30 mm

**03 Back systems**

CODE	DESCRIPTION
737	Standard back, back support height 400 mm
736	Adjustable back, back support height 400 mm
26	30° Back adjustment through pneumatic spring  for back system CODE 737 and 736

**04 Arm supports and arm pads**

CODE	DESCRIPTION
106	Standard arm support with upholstered arm pad and black clothes guard

☐ Please select by ticking the relevant boxes☒ Special construction/individual solution o.r. - on request☐ Standard version: If no other selection is made, the standard version will be delivered.**05 Leg supports and footrests**

CODE	DESCRIPTION
93/798	Legrest, lift off and swing away, divided footrest, hight adjustable, flip up, angle adjustable
92/798	Legrest, swing-up, lift off and swing away, divided footrests, hight adjustable, flip up and angle adjustable

**06 Steering wheel tyres**

CODE	DESCRIPTION
157/481	Pneumatic tyres 200 x 50 mm (8")
157/900	Puncture safe tyres 200 x 50 mm (8")

**07 Tyres drive wheels**

CODE	DESCRIPTION
482/815	Pneumatic tyres 20 x 60 mm (12.5")
815/904	Puncture safe tyres 20 x 60 mm (12.5")

**08 Batteries and Chargers**

CODE	DESCRIPTION
4696	Maintenance-free batteries 38Ah (5h) / 45 Ah (20h)
731	Battery charger 24 V / 6 A

**09 Operating modules and brackets**

CODE	DESCRIPTION
408	Operating module VR 2
60	Mounted left
74	Assembly right
4859	Operating module bracket, length adjustable
4860	Operating module bracket with parallelogram length and height adjustable

☐ Please select by ticking the relevant boxes☒ Special construction/individual solution o.r. - on request☐ Standard version: If no other selection is made, the standard version will be delivered.**10 Lighting**

CODE	DESCRIPTION
932	Active lighting

**11 Frame- and cover colours**

CODE	DESCRIPTION
4919	Silverline
287	Seat and back strap black

**12 Options around the wheelchair**

CODE	DESCRIPTION
970	Crutch holder
833	Safety belt

**13 software**

CODE	DESCRIPTION
985	Series programming (standard)