

iTRAVEL

1.054 Stock version



Allowed user weight

120 kg





Adjustable backrest



Comfort seat system



Flip-up armrests



Small pack size

Travel Wheelchair

ADAPTATION

- Weighing in at 22 kg, this is the lightweight among the electric wheelchairs
- The easy-fold system makes the iTravel the perfect travelling companion
- A walking stick holder is mounted as standard

OPERATION

- Very easy to use

ORDER

COST ESTIMATE

CUSTOM DESIGN

Please send the order form by e-mail to:

info@meyragroup.com

Customer service: +49 571 93292 555

Billing address / Customer ref. no.:

Delivery address:

Tax-ID

Consignment / Comments:

ORTOPEDIA

EINE MARKE VON MEYRA



Basic equipment

- Swing-up armrests
- Easy operation with Speed Preselection
- Battery charger

- Sensor
- Seat foldable
- Battery removable with additional charging socket



Allowed user weight

120 kg



Hints

Our Terms and Conditions as stated at www.meyra.de/AGB shall apply

Valid from 01.01.2025 – 31.12.2025.

Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de

You avoid queries and incorrect orders thanks to the buildability check



Technical specifications

Length in mm	970
Width in mm	560
Seat depth in mm	400
Seat widths in mm	440
Height in mm	870
Back support height in mm	460
Height when folded in mm	300
Footplate to seat (lower leg length) in mm	440
Empty weight from	22 kg
User weight	120 kg
Driving speed	6 km/h
Castor wheels diameter	200
Drive wheels diameter	229
Ground clearance in mm	70
Max. obstacle height in mm	40
Turning radius in mm	970
Continuous motor performance	200 W
Peak motor performance	800 W
Batteries	10,4 Ah with lithium battery
Battery charger	2 A
Electronic	50 A
Empty weight with lithium batteries	23 kg
Permissible overall weight	150 kg

Please enter your quantities in the boxes

01 iTravel E-Mobile

Item	DESCRIPTION	Chassis colour
9105414	MEYRA iTRAVEL Lithium Ionen Battery	White

02 Accessories

Item	DESCRIPTION
10542837	Additional lithium battery <i>*Nur i.V. mit Modell 1.054 möglich</i>

USER DATA

First name

Name

Street

ZIP CODE

City

Country

Birthday

Diagnosis 1

Diagnosis 2

Tel.

Mobil

E-mail

Order number

Health insurance fund / cost unit

Specialist trade Contact person/phone.

Supply specialist MEYRA

Approval info from week

ANATOMICAL DATA (IN MM)

1) Seat depth

2) Seat width

3) Chest width

4) Chest depth

5) Elbow height

6) Armpit height

7) Back height


8) Head height

9) Forearm length

10) LLL

Body height

Weight (kg)



Remarks

I accept the Terms and Conditions of Meyra Distribution GmbH (valid as of January 2023)

<https://www.meyra.com/gtc/>

Date

City

Signature