

MICRO ELECTRO MOBIL

1.064 Stock version

 **Allowed user weight**
120 kg





500 mm turning radius



Dismantles into three parts



Simple operation



Mid-wheel drive

A real space saver

ADAPTATION

- Compact dimensions
- Lightweight construction, can be splitted in three parts, thus allowing uncomplicated transport

OPERATION

- Very easy to use
- Easy turning through the mid-wheel drive and a turning radius of 500 mm

ORDER

COST ESTIMATE

CUSTOM DESIGN

Please send the order form by e-mail to:

info@meyragroup.com

Customer service: +49 571 93292 555

Billing address / Customer ref. no.:

Delivery address:

Tax-ID

Consignment / Comments:

MEYRA



Basic equipment

- Swing-up armrests
- self-explanatory Display
- Battery charger
- Sensor

- Seat rotatable
- Seat height-adjustable
- Seat detachable
- Battery box removable with additional charging socket



Allowed user weight

120 kg



Hints

Our Terms and Conditions as stated at www.meyra.de/AGB shall apply

Valid from 01.01.2025 – 31.12.2025.

Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de

You avoid queries and incorrect orders thanks to the buildability check

**Technical specifications**

Length in mm	760
Width in mm	580
Seat height adjustable in mm	370 - 450
Seat depth in mm	400
Seat width adjustable in mm	470 - 570
Back support height in mm	400
Transport height in mm	470
Castor wheels diameter	110
Drive wheels diameter	220
Ground clearance in mm	50
Max. obstacle height in mm	15
Turning radius in mm	500
Continuous motor performance	200 W
Peak motor performance	800 W
Batteries	20,8 Ah lithium battery / 14 Ah lead battery
Battery charger	2 A
Electronic	50 A
Range with lithium batteries	17 km
Range with lead batteries	10 km
Permitted inclination/slopes	3° (6°)
Empty weight with lithium batteries	42 kg
Empty weight with lead batteries	47 kg
User weight	120 kg
Permissible overall weight	170 kg
Driving speed	6 km/h

Please enter your quantities in the boxes

01 1.064 E-Mobile

Item	DESCRIPTION	Chassis colour
9106404	ELEKTROMOBIL lead acid batteries	White
9106414	ELEKTROMOBIL lithium ion batteries	White
9106424	ELEKTROMOBIL up lead acid batteries	White
9106434	ELEKTROMOBIL up lithium ion batteries	White

02 Options

Item	DESCRIPTION
1083217	Shopping basket with holder
1081860	Soft wheels front twice and back once

USER DATA

First name

Name

Street

ZIP CODE

City

Country

Birthday

Diagnosis 1

Diagnosis 2

Tel.

Mobil

E-mail

Order number

Health insurance fund / cost unit

Specialist trade Contact person/phone.

Supply specialist MEYRA

Approval info from week

ANATOMICAL DATA (IN MM)

1) Seat depth

2) Seat width

3) Chest width

4) Chest depth

5) Elbow height

6) Armpit height

7) Back height


8) Head height

9) Forearm length

10) LLL

Body height

Weight (kg)



Remarks

I accept the Terms and Conditions of Meyra Distribution GmbH (valid as of January 2023)

<https://www.meyra.com/gtc/>

Date

City

Signature